

## *Certificate of Analysis*

Product Name: U.C.R.P. - Orig:USA  
Product Numbers: 100301  
Lot Number: 281760

Production Date: 02-MAR-13  
Expiration Date: 2015-07-31

| Specification            | Specification Limit                                                                                    | Result   |
|--------------------------|--------------------------------------------------------------------------------------------------------|----------|
| Appearance Cake          | Yellow to buff colored solid                                                                           | Complies |
| Appearance Reconstituted | Yellow to straw colored, clear to hazy,<br>free of particulates when viewed against<br>dark background | Complies |
| Solubility               | Less than 2 minutes                                                                                    | Complies |
| Moisture                 | $\leq 3$ %                                                                                             | 1        |

For in-vitro diagnostic use.

See directional insert for assay values.

Each unit of source material used in the preparation of this product has been tested by an FDA licensed method and found nonreactive for HBsAg and negative for antibodies to HIV and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS, or other infectious diseases. This product, like all materials of human origin, should be handled as potentially infectious biological material.



Quality Representative

Date: 2013-03-05