



Certificate of Analysis

Globe Scientific certifies that this product, made from resin material, has been manufactured under controlled conditions and is in compliance with all applicable procedures and specifications.

Item #: 1334Y

Lot #: 102493-19-011816

Expiration Date: N/A

Product Descriptions Microscope slide, yellow frost beveled edge, clip corner, 72/bx

Component Material: Soda Lime Glass

Manufacturing Standard(s): ISO 9001:2008, ISO 13485:2003, CE 93/42/EEC, CE (IVD) 98/79/EEC

Globe Scientific certifies that this product has been tested and is completely free of the following where applicable:

Statements:

Transportation and storage conditions considerably affect the quality and performance of microscope slides. Microscope slides should be stored in a dry place at constant room temperature. Slides should be stored on a pallet or shelf and not directly on a floor. Do not store slides and solvents together as the slide surface can turn hydrophobic. If the temperature is constantly warming and cooling, condensation may occur, which can cause sticking of the slides. We also recommend, when shipping to a customer, instruct them to allow the slides to reach room temperature prior to opening the cartons and using them.

Since soda lime glass is subject to an aging process, we further recommend continuously rotating stock so the oldest slides are always used first. This will ensure that all slides being used are of the highest quality.

Standard shelf life for soda lime glass microscope slides is approximately 2 years if the above transportation and storage conditions are followed.

A blue ink signature of David Ackley, written in a cursive style.

David Ackley
Vice President, Product Development

07/25/18

Date

