



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications

Microorganism Name: Clostridium sporogenes
Catalog Number: 0487
Lot Number: 487-263
Reference Number: ATCC® 11437™*
Purity: Pure
Passage from Reference: 3

Expiration Date: 2018/12/31

Release Information:

Quality Control Technologist: Mary L Bowman

Release Date: 2017/2/21

Performance

Macroscopic Features:

Large, irregularly circular, with coarse rhizoid edge, have a raised yellowish-gray center and flattened periphery, beta hemolytic, adheres to agar.

Medium:

A/R SBAP

Microscopic Features:

Gram positive straight rods. Spores, when seen, are oval, central or subterminal and distend the cell slightly.

Method:

Gram Stain (1)

ID System: MALDI-TOF

See attached ID System results document.

Other Features/ Challenges: Results

Egg yolk agar: lecithinase - negative

Egg yolk agar: lipase - positive

Amanda Kuperus

Quality Control Manager

AUTHORIZED SIGNATURE

Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.



(1) These tests are accredited to ISO/IEC 17025:2005.

TESTING CERT #2655.01



Meaning of Score Values

Range	Interpretation	Symbols	Color
2.00 - 3.00	High-confidence identification	(+++)	green
1.70 - 1.99	Low-confidence identification	(+)	yellow
0.00 - 1.69	No Organism Identification Possible	(-)	red

Meaning of Consistency Categories (A - C)

Category	Interpretation
(A)	High consistency: The best match is a high-confidence identification. The second-best match is (1) a high-confidence identification in which the species is identical to the best match, (2) a low-confidence identification in which the species or genus is identical to the best match, or (3) a non-identification.
(B)	Low consistency: The requirements for high consistency are not met. The best match is a high- or low-confidence identification. The second-best match is (1) a high- or low-confidence identification in which the genus is identical to the best match or (2) a non-identification.
(C)	No consistency: The requirements for high or low consistency are not met.

Sample Name: Clostridium sporogenes
 Sample Description: 0487
 Sample ID: 487-263
 Sample Creation Date/Time: 2017-02-21T08:38:05.294 MB
 Applied MSP Library(ies): BDAL, Mycobacteria Library (bead method), Filamentous Fungi Library 1.0, Listeria

Sample Name	Sample ID	Organism (best match)	Score Value
C3 (+++)(A)	487-263	Clostridium sporogenes	2.07

Comments:

Clostridium sporogenes is closely related and shows very similar spectra to the strains of the highly pathogenic Clostridium botulinum groups A, B and F. Clostridium botulinum is not included in the MALDI Biotyper database. For differentiation an adequate identification method has to be selected by an experienced professional.