



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Eikenella corrodens Catalog Number: 0189 Lot Number: 189-35** Reference Number: ATCC® BAA-1152™* Purity: Pure Passage from Reference: 3	Expiration Date: 2019/12/31 Release Information: Quality Control Technologist: Carol J Stanoch Release Date: 2018/2/9
---	--

Performance	
Macroscopic Features: Medium, circular, convex, entire edge, light tan, glistening, "bleach-like" odor; as culture ages, margins flatten out.	Medium: Chocolate
Microscopic Features: Gram negative straight rod with rounded ends.	Method: Gram Stain (1)
ID System: Vitek NH (1) See attached ID System results document.	Other Features/ Challenges: Results (1) Oxidase (Kovacs): positive (1) Catalase (3% Hydrogen Peroxide): negative MacConkey Agar (AER/35 C/48 Hr): no growth  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE

**Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.



bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Feb 9, 2018 07:36 CST
Printed by: cjs

Isolate Group: 189 35-1

Bench: CS

Card Type: NH Testing Instrument: 000017919CE3 (15220)

Bionumber: 0622002400

Organism Quantity:

Comments:	

Identification Information		Card: NH	Lot Number: 2450361203	Expires: Nov 22, 2018 12:00 CST
		Completed: Feb 8, 2018 14:58 CST	Status: Final	Analysis Time: 6.00 hours
Selected Organism		97% Probability Eikenella corrodens		
		Bionumber: 0622002400	Confidence: Excellent identification	
SRF Organism				
Analysis Organisms and Tests to Separate:				
Analysis Messages:				
Contraindicating Typical Biopattern(s)				
Eikenella corrodens dMLT(16),				

Biochemical Details

1	ArgA	-	2	GGT	-	3	LysA	-	4	dGAL	-	5	LeuA	+	6	ELLM	+
7	PheA	-	8	ProA	+	10	PyrA	-	13	TyrA	-	15	APPA	+	18	dGLU	-
19	GLYG	-	20	dMNE	-	22	dMAL	-	28	SAC	-	33	NAG	-	36	URE	-
39	BGALi	-	40	ODC	+	41	AARA	-	45	PVATE	-	46	PHC	-	47	dMLT	+
51	MTE	-	52	IGLM	(-)	59	PHOS	-	61	dRIB2	-	62	OPS	-	64	dXYL	-

Installed VITEK 2 Systems Version: 07.01

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: