



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Clostridium perfringens Catalog Number: 0318 Lot Number: 318-154 Reference Number: ATCC® 13124™* Purity: Pure Passage from Reference: 2	Expiration Date: 2019/1/31 Release Information: Quality Control Technologist: Tracy A Blenker Release Date: 2017/3/28
---	--

Performance	
Macroscopic Features: Gray, opaque, low, flat, somewhat rhizoid: tend to spread but not swarm, double zone of hemolysis	Medium: A/R SBAP
Microscopic Features: Gram positive straight rods. Spores are rarely seen; if present they are large oval, central or subterminal and distend the cell.	Method: Gram Stain (1)
ID System: Vitek ANC (1) See attached ID System results document.	Other Features/ Challenges: Results Egg yolk agar: lecithinase - positive Egg yolk agar: lipase - negative Reverse CAMP reaction: positive  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE

Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.



(1) These tests are accredited to ISO/IEC 17025:2005.

TESTING CERT #2655.01

bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Mar 20, 2017 08:55 CDT
Printed by: tab

Isolate Group: 318 154-1

Bench: TB

Card Type: ANC Testing Instrument: 000017919CE3 (15220)

Bionumber: 0427505250511

Organism Quantity:

Comments:	

Identification Information	Card: ANC	Lot Number: 2440202103	Expires: Jun 16, 2018 13:00 CDT
	Completed: Mar 17, 2017 14:55 CDT	Status: Final	Analysis Time: 6.00 hours
Selected Organism	98% Probability Clostridium perfringens Bionumber: 0427505250511 Confidence: Excellent identification		
SRF Organism			
Analysis Organisms and Tests to Separate:			
Analysis Messages:			
Contraindicating Typical Biopattern(s)			

Biochemical Details

4	dGAL	-	5	LeuA	-	6	ELLM	-	7	PheA	(-)	8	ProA	-	10	PyrA	+
11	dCEL	-	13	TyrA	+	15	APPA	-	18	dGLU	+	20	dMNE	+	22	dMAL	+
28	SAC	+	30	ARB	-	33	NAG	+	34	BGLUi	-	36	URE	-	37	BGURi	-
39	BGALi	+	41	AARA	-	42	AGALi	+	43	BMAN	-	44	ARG	+	45	PVATE	-
51	MTE	+	53	ESC	-	54	BdFUC	+	55	BNAGi	-	56	AMANi	-	57	AIFUC	-
59	PHOS	+	60	IARA	-	61	dRIB2	+	62	OPS	+	63	AARAF	-	64	dXYL	-
	GRAM	+		MORPH	-		AERO	-									

Installed VITEK 2 Systems Version: 07.01

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: