



DIESSE DIAGNOSTICA SENESE SpA
Via delle Rose 10, 53035 Monteriggioni (SI), Italy
Tel. ++39-0577-587111 fax ++39-0577-318690

Sede Legale e Amministrazione
Via A. Solari 19 20144 Milano Italy
Tel. ++39-02-4859121 r.a. fax ++39-02-48008530

CERTIFICATE OF ANALYSIS

CODE	PRODUCT NAME	LOT	EXPIRY DATE	PREP DATE
10435	ESR Control 4 x 9 ml	239	2017-09	2016-09
10436	ESR Control 2 x 9 ml	239	2017-09	2016-09

Level I (NORMAL Level)	Test conditions	VES 1st hour	Range of acceptability
VES-MATIC CUBE 30	Temperature OFF	6	1 - 11
VES-MATIC CUBE 80 - 200	Temperature OFF	7	2 - 12

Level II (ABNORMAL Level)	Test conditions	VES 1st hour	Range of acceptability
VES-MATIC CUBE 30	Temperature OFF	73	48 - 98
VES-MATIC CUBE 80 - 200	Temperature OFF	50	25 - 75

Date: 19/12/2016


Alessandra Cianchi
Signature (Head Quality Control)