

# CERTIFICATE OF CONFORMITY

## Certificate for Lot Number: 52748

Product Code: TS/72-MX      Product: Mixed Protect, Bulk Bag - Sale size 1000 (Pack Size 500)

Date Of Manufacture: 04/2017    Expiry: 2019-04      Shelf Life: 24 months

This is to certify that the product detailed on this certificate conforms to Technical Service Consultants Ltd procedures and specifications and has performed the tests detailed below within documented parameters.

### Sterilisation Reference, Isotron PLC Number

S11845389-2-1 - 5178822 : Dose 23.1 - 31.5 kGy, S11850088-5-1 - 5206482 : Dose 27.2 - 37.2 kGy, S11853980-6-1 - 5240133 : Dose 25.7 - 35.1 kGy, S11859151-6-1 - 5256902 : Dose 23.7 - 32.4 kGy

For Irradiated products batch release is controlled by the verification of the radiation dose throughout the load.

### Autoclave Reference, Time Manufactured

For autoclaved products, either a temperature datalogger or visual sterilisation indicator are used to confirm that sterility was achieved.

### Physical and Microbiological Assurance

| Tests Performed                                                                  | Results Obtained |
|----------------------------------------------------------------------------------|------------------|
| Sterility                                                                        | Pass             |
| Product Specification                                                            | Pass             |
| Recovery of <i>S. pyogenes</i> ATCC 19615 on Chocolate Agar/CO <sub>2</sub> /36± | Pass             |
| Recovery of <i>H. Influenzae</i> ATCC 49247 on Chocolate Agar/CO <sub>2</sub> /3 | Pass             |
|                                                                                  |                  |

**Overall Result: Pass**

Certificate Issued by: Joanne Walsh    Certificate authorised on 28.04.2017

Electronically generated on 16th October 2017 - valid without signature

Sent to Vickie (vickie.bruckshaw@tscswabs.co.uk)