



CERTIFICATE OF ANALYSIS

KIT NAME: Mythic 18-22 FLUSH-CLEANER

CAT. NO: HM18-011-025

LOT NO: 708-01

EXPIRY DATE: 2018-09

DATE OF QC: 25.09.2017

<i>PARAMETER</i>	<i>NORM</i>	<i>RESULTS</i>
Appearance	clear	clear
Colour	yellow	yellow
Odour	specific smell (chlorine)	specific smell (chlorine)
Active chlorine concentration (g/dm ³)	32,5 (30,0-35,0)	34,67

We hereby certify that this lot has been tested and was found to meet our specifications.
The production is in compliance to Quality Management System ISO 13485:2003.

Sign.: Anna Spryszak
Quality Assurance Department

This certificate is a computer printout and has therefore not been signed by hand.

PRINTED VERSION IS NOT SUPERVISED, CHECK TOPICALITY OF DOCUMENT BEFORE USE

Document number	Document version	Date of revision	Prepared by	Approved by
FOR039	2	04.11.2015	RW	JP