

CERTIFICATE OF ANALYSIS

Reagent: Microalbumin
Catalogue Number: 252-20
Lot Number: 50441
Expiry Date: 2018-06-30
Storage Conditions: 2-8°C

A. PHYSICAL CHARACTERISTICS

PARAMETER	SPECIFICATION	RESULT
Optical Density of Reagent at 340 nm at 18-26°C: - Buffer Reagent (R1) - Antibody Reagent (R2)	≤ 0.030 0.100 - 0.383	0.008 0.175
pH of Reagent at 25°C - Buffer Reagent (R1) - Antibody Reagent (R2)	7.0 - 8.0 7.0 - 8.0	7.7 7.7
Linearity	30 mg/dL	Pass
Performance	Values obtained for normal and abnormal level control sera are within the ranges established for the method.	

B. CHEMICAL CHARACTERISTICS

Reagent Component

Major Active Ingredients

Microalbumin Buffer Reagent (R1)
Preservative

tris (hydroxymethyl) aminomethane 100 mmol/L
sodium azide

Microalbumin Antibody Reagent (R2)
Preservative

tris (hydroxymethyl) aminomethane 100 mmol/L
anti-human albumin, in goat anti serum 20%
sodium azide

Intended Use

For the quantitative determination of low levels of albumin in urine. For IN VITRO diagnostic use.

The above reagent has passed all of our quality control procedures and has been released for distribution.

21-Sep-2017
Date

S. M. J. 00
Quality Control

Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Tel: 800-332-1042 Fax: 800-762-6311

Sekisui Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
Kent, ME16 0LQ UK
Tel: +44 (0)1622 607800 Fax: +44 (0)1622 607801

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, Prince Edward Island
C1E 2B9 Canada
Tel: 800-565-0265 Fax: 902-628-6504