

## *Certificate of Analysis*

Product Name: Calcium Chloride, 0.02M - Orig:USA  
Product Number: 100304  
Lot Number: 218365


Expiration Date: 09/30/2020  
Storage: 2 - 8°C

| Specification                           | Specification Limit                                                 | Result   |
|-----------------------------------------|---------------------------------------------------------------------|----------|
| Appearance                              | Clear, colorless solution essentially<br>free of particulate matter | Complies |
| Performance v. Previous Lot: Difference | $\leq 3$ %                                                          | 0        |

For in-vitro diagnostic use.

Reactive Component: Calcium Chloride

WARNING: Calcium chloride solution contains sodium azide. Sodium azide under acidic conditions yields hydrazoic acid, an extremely toxic compound. Dilute with running water before discarding, and then flush with large volumes of water. These precautions are recommended to avoid deposits in metal piping in which explosive conditions may develop.



Quality Representative

Date: 2017-10-25