



# CERTIFICATE OF CONFORMITY

## Certificate for Lot Number: 43424

Product Code: TS/17-B      Product: Polypropylene/Viscose Swab in Peel Pouch (Pack Size 1000)

Date Of Manufacture: 12/2013    Expiry: 2015-12      Shelf Life: 24 months

This is to certify that the product detailed on this certificate conforms to Technical Service Consultants Ltd procedures and specifications and has performed the tests detailed below within documented parameters.

Sterilisation Reference, Isotron PLC Number

S11086716-5-1 - 663584 : Dose 32.1 - 40.7 kGy

For Irradiated products batch release is controlled by the verification of the radiation dose throughout the load.

Autoclave Reference, Time Manufactured

For autoclaved products, either a temperature datalogger or visual sterilisation indicator are used to confirm that sterility was achieved.

### Physical and Microbiological Assurance

| Tests Performed | Results Obtained |
|-----------------|------------------|
| Sterility       | Pass             |
|                 |                  |
|                 |                  |
|                 |                  |
|                 |                  |

**Overall Result: Pass**

Certificate Issued by: Joanne Walsh    Certificate authorised on 16.12.2013

Electronically generated on 3rd January 2014 - valid without signature

Sent to vickie (vickie.bruckshaw@tscswabs.co.uk)