

Certificate of Analysis

Product Description: **ACCURUN 190 Series 1000 Anti-Trypanosoma cruzi (Chagas) Positive Control**
 Product Code: **2015-0098** Lot/Batch Number: **10219885**
 Storage: **2-8°C** Expiration Date: **20 JUL 2018**

Test Results: The Abbott ARCHITECT Chagas test method was used to evaluate ACCURUN 190 Series 1000 Anti-Trypanosoma cruzi (Chagas) Positive Control. The test values in the table below represent the mean of 4 replicates.


			Acceptance Range (s/co)	
Description	Batch	Test Value (s/co)	Lower Limit	Upper Limit
ACCURUN 190 Series 1000	10219885	2.1	1.5	2.9

Raw Material: Human Plasma reactive for Anti-Trypanosoma cruzi (Chagas).

Matrix: Defibrinated human plasma negative for HBsAg and antibodies to HIV-1 and HIV-2, HTLV I/II, and HCV.

Preservatives: 0.1% ProClin300


Approval:



Prepared By

5 APR 17

Date



QA Verified By

05 Apr 17

Date