

ACON BIOTECH (HANGZHOU) CO.,LTD

CERTIFICATE OF ANALYSIS

PRODUCT NAME: Urinalysis Reagent Strips

PURCHASE NO.: 43040001

BATCH NO.: URS6100108

EXPIRATION DATE: 2018-08

QUANTITY: 100x100T

ANALYSIS		SPECIFICATIONS (Color card)	TEST RESULT
Calcium	MAS urine control level 1	≥ 4 mg/dL <10mg/dL	Complies
	40mg/dL	>30mg/dL ≤ 40 mg/dL	Complies
Blood	Neg. urine	Negative	Complies
	Yihua Urine controls	>+ <+++	Complies
pH	5.0	≥ 5.0 <6.0	Complies
	9.0	>8.0 ≤ 9.0	Complies
Creatinine	MAS urine control level 1	≥ 10 mg/dL <50mg/dL	Complies
	MAS urine control level 2	>100mg/dL ≤ 300 mg/dL	Complies
Nitrite	Neg. urine	Negative	Complies
	Yihua Urine controls	+	Complies
Ketone	Neg. urine	Negative	Complies
	Yihua Urine controls	>+ <+++	Complies
S.G	1.000	≥ 1.000 <1.005	Complies
	1.025	>1.020 <1.030	Complies
Ascorbic acid	Neg. urine	Negative	Complies
	20mg/dL	>+ <+++	Complies
Glucose	Neg. urine	Negative	Complies
	Yihua Urine controls	>+ <++++	Complies
Bilirubin	Neg. urine	Negative	Complies
	Yihua Urine controls	>+ $\leq$ +++	Complies
Protein	Neg. urine	Negative	Complies
	Yihua Urine controls	>+ <+++	Complies
Albumin	MAS urine control level 1	≥ 10 mg/L <30mg/L	Complies
	MAS urine control level 2	>80 mg/L ≤ 150 mg/L	Complies
Urobilinogen	Neg. urine	Negative	Complies
	4~8mg/dL	>2mg/dL <12mg/dL	Complies
Leukocyte	Neg. urine	Negative	Complies
	Yihua Urine controls	>+ <+++	Complies

QC Supervisor: Ivy Fang
Signature: Ivy Fang
Date: Oct 26, 2016