

DIESSE DIAGNOSTICA SENESE SpA

53035 Monteriggioni (SI) – Italia – Via delle Rose 10
Tel +39-0577-587111 Fax +39-0577-318690

**Sede Legale e Amministrazione**

20144 Milano (Italy) Via A. Solari 19 sc. 6
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

CERTIFICATO DI CONTROLLO QUALITÀ/QUALITY CONTROL CERTIFICATE**PRODOTTO/PRODUCT**

ESR CONTROL CUBE

REF:

10435 (4x9 ml)

10436 (2x9 ml)

REAGENTI/REAGENTS

	Kit	Level 1	Level 2
Lot:	311	311	311
	2017-03	2017-03	2017-03
	2018-03	2018-03	2018-03

RISULTATI/RESULTS

	Level 1		Level 2		Temp.
	Risultato/ Result	Range	Risultato/ Result	Range	
VES-Matic Cube 30	7	2-12	83	58-108	OFF
VES-Matic Cube 80/VES-Matic Cube 200	6	1-11	61	36-86	OFF

Data/Date

11/04/2017

Testato da/Tested by:
Katia Tanzini**Approvato da/Approved by:**
Alessandra Clanchi (Responsabile CQ/QC Head)