

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI) – Italia – Via delle Rose 10
Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

Pag 1/1

QUALITY CONTROL CERTIFICATE**KIT NAME: snRNP-C - CHORUS**

Code 86092	
Lot No. 302	Expiry 2018-01

REAGENTS

Label	Description	Lot.	Exp	No
D7441	snRNP-C	302	2018-01	6X6
B7441	CALIBRATORE snRNP-C	302	2018-01	1
C7441	CONTROL SERUM snRNP-C	302	2018-01	1

CONTROLS


	Found results	Expected results
Calibrator (mOD)	420	≥ 100
Control Serum (AU/mL)	53,2	27,0-58,0

CALIBRATOR

1	9	1	3	0	2	6	3	2	0	3	5	3	8	5	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader


YES

Operator Signature *Antonella Leo***CONTROL SERUM**

1	9	1	3	0	2	6	1	8	3	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES

Operator Signature *Antonella Leo*

The test has been carried out in accordance with the package insert

TESTED BY: Antonella Leo**APPROVED** RESPONSIBLE HEAD, Quality Control Reagents *Antonella Leo***DATE** 30/01/2017**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.