

Certificate of Analysis

Product Name: Factor VIII Deficient Plasma - Orig:USA
Product Number: 100800
Lot Number: 823109

Expiration Date: 04/30/2019
Storage: 2 - 8°C

Specification	Specification Limit	Result
Appearance Cake	Yellow to buff colored solid	Complies
Appearance Reconstituted	Clear to slightly cloudy, straw colored liquid free of particulates.	Complies
Factor VIII Activity	< 1.0 %	0.8
APTT	>= 60 seconds	Complies
Prothrombin Time	10 - 13 S	13
Moisture	<= 3 %	1
Factor VIII	No Significant Correction*	84 S
Factor IX	Significant Correction*	42 S
Factor XI	Significant Correction*	45 S
Factor XII	Significant Correction*	40 S

For in-vitro diagnostic use.

* Correction when mixed 1:1 with product noted.

All individual donor units comprising this lot of Factor VIII Deficient Plasma have been tested and found nonreactive for HBsAg and HCV and negative for antibodies to HIV. As with all blood products, this product should be handled as potentially infectious biological material.



Quality Representative

Date: 2016-04-19