

Certificate of Analysis

Product: SeraCon II**Product Code: 22000****Lot Number: 10E0608****Storage: -20°C****Manufacture Date: 24SEP10**

<u>Analyte</u>	<u>Results</u>	<u>Assay Lower Limit</u>	
Albumin	3.8	g/dL	
Calcium	5.9	mg/dL	
Cholesterol	135	mg/dL	
Triglyceride	79	mg/dL	
** BUN	0	mg/dL	2 mg/dL
** Creatinine	0.1	mg/dL	0.2 mg/dL
Glucose	8	mg/dL	
T. Protein	6.2	g/dL	
** Uric Acid	0.3	mg/dL	1.5 mg/dL
Sodium	100	mmol/L	
Potassium	4.6	mmol/L	
Chloride	79	mEq/L	
** Phosphorus	0.7	mg/dL	1 mg/dL
Bacteria	1	cfu/mL	
pH	7.4		

** Reported test result is below lower limit of assay linearity

This product was manufactured from human plasma and each donor unit used was tested by an FDA approved method for the presence of antibody to HIV 1/2, antibody to HCV, HIV-Ag and/or HIV NAT as well as for HBsAg, and found to be negative.

Prepared By: Charlene Lachapelle
Charlene Lachapelle, QA Specialist II

01 Oct 10
Date

QA Review By: Kerz Poku
Kerz Poku, QA Specialist II

04 Oct 10
Date