

RANDOX

CERTIFICATE OF ANALYSIS

This is to certify that the following kit was tested by a
Randox Testing Laboratory and approved by the QC Department.

DATE OF MANUFACTURE: November 2016
KIT DESCRIPTION: IA PREMIUM PLUS TRI-LEVEL
CATALOGUE NO: IA3112
LOT NUMBERS: 398901 (1619EC, 1644EC, 1646EC)
EXPIRY DATE: 2019-04-28
SYSTEM OF ANALYSIS: Siemens Immulite 2000 Xpi

TEST CONTROLS OR VIALS

ANALYTE	LOT NUMBER	RANGE	TARGET	RESULT	PRECISION (if applicable)
CEA (ug/l)	1619EC	3.82-5.72	4.77	4.8	
				4.8	
	1644EC	19.1-28.7	23.9	23.3	
				23.2	
	1646EC	59.8-89.6	74.7	76.6	
				80.8	
C-peptide (pmol/l)	1619EC	287-478	382	387	
				394	
	1644EC	1070-1780	1420	1427	
				1446	
	1646EC	1330-2210	1770	1784	
				1821	
DHEAS (umol/l)	1619EC	1.64-2.46	2.05	1.9	
				1.9	
	1644EC	9.6-14.4	12	11.8	
				12.2	
	1646EC	16.4-24.6	20.5	21.5	
				22.1	

QC TESTED BY: Davina Graham
DATE: 07 Nov 2016
RESULT (PASS/FAIL): pass
QC APPROVED BY: Adele Martin.

Randox Laboratories Limited, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom
T +44 (0) 28 9442 2413 F +44 (0) 28 9445 2912 E marketing@randox.com I www.randox.com

Directors: Dr S P Fitzgerald BSc PhD DSc (Hons) FIAE Man. R J Fitzgerald, W Fitzgerald OBE, R P Kelly BSc FCA.
Randox Laboratories Limited is a company registered in Northern Ireland. Company Number NI15738. VAT Number GB 353 0304 00.

FORM NO. 6521-QC
REVISION (4)
14 April 2016