



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications

Microorganism Name: Escherichia coli
Catalog Number: 0335
Lot Number: 335-389**
Reference Number: ATCC® 25922™*
Purity: Pure
Passage from Reference: 3

Expiration Date: 2020/7/31

Release Information:

Quality Control Technologist: Christine Condon

Release Date: 2018/8/16

Performance

Macroscopic Features:

2 colony types, both are gray & beta hemolytic; one is circular to irregular, convex, slightly erose edge & smooth; other is larger, irregular, low convex, erose edge & rough

Medium:

SBAP

Microscopic Features:

Gram negative straight rod

Method:

Gram Stain (1)

ID System: Vitek GN (1)

See attached ID System results document.

Other Features/ Challenges: Results

(1) Oxidase (Kovacs): negative

Beta-glucuronidase (E. coli Broth w/MUG): positive

(1) Ampicillin (10 mcg - Disk Susceptibility): 16 - 22 mm

(1) Gentamicin (10 mcg - Disk Susceptibility): 19 - 26 mm

(1) SXT (1.25/23.75 mcg - Disk Susceptibility): 23 - 29 mm

Amanda Kuperus

Quality Control Manager

AUTHORIZED SIGNATURE

****Disclaimer:** The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.



bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Aug 2, 2018 17:27 CDT
Printed by: cmc
Report Version: 2 of 2

Isolate Group: 335 389-1

Bench: CC

Last Updated: Aug 2, 2018 17:27 CDT

By: cmc

Card Type: GN Testing Instrument: 000017919CE3 (15220)

Bionumber: 0405611560567601

Organism Quantity:

Comments:	
	QC OK

Identification Information		Card: GN	Lot Number: 2410576403	Expires: Jun 25, 2019 13:00 CDT
		Completed: Aug 1, 2018 17:46 CDT	Status: Final	Analysis Time: 4.75 hours
Selected Organism		98% Probability Escherichia coli		
		Bionumber: 0405611560567601	Confidence: Excellent identification	
SRF Organism				
Analysis Organisms and Tests to Separate:				
Analysis Messages:				
Contraindicating Typical Biopattern(s)				
Escherichia coli dTAG(22),				

Biochemical Details																	
2	APPA	-	3	ADO	-	4	PyrA	-	5	IARL	-	7	dCEL	-	9	BGAL	+
10	H2S	-	11	BNAG	-	12	AGLTp	-	13	dGLU	+	14	GGT	-	15	OFF	+
17	BGLU	-	18	dMAL	+	19	dMAN	+	20	dMNE	+	21	BXYL	-	22	BAlap	-
23	ProA	+	26	LIP	-	27	PLE	-	29	TyrA	+	31	URE	-	32	dSOR	+
33	SAC	-	34	dTAG	+	35	dTRE	+	36	CIT	-	37	MNT	-	39	5KG	-
40	ILATk	+	41	AGLU	-	42	SUCT	+	43	NAGA	-	44	AGAL	+	45	PHOS	+
46	GlyA	+	47	ODC	+	48	LDC	+	53	IHISa	-	56	CMT	+	57	BGUR	+
58	O129R	-	59	GGAA	-	61	IMLTa	-	62	ELLM	+	64	ILATa	-			

Installed VITEK 2 Systems Version: 07.01

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: