



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Acinetobacter baumannii Catalog Number: 0119 Lot Number: 119-32** Reference Number: ATCC® BAA-747™* Purity: Pure Passage from Reference: 3	Expiration Date: 2019/12/31 Release Information: Quality Control Technologist: Tracy A Blenker Release Date: 2018/2/13
--	---

Performance	
Macroscopic Features: Large, circular, convex, entire edge, cream, glistening Microscopic Features: Short, fat gram negative rods, singly and in pairs	Medium: SBAP Method: Gram Stain (1)
ID System: Vitek GN (1) See attached ID System results document.	Other Features/ Challenges: Results (1) Oxidase (Kovacs): negative Growth at 41 C: positive Growth at 44 C: positive  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE

**Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.



bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Feb 12, 2018 16:44 CST
Printed by: tab

Isolate Group: 119 32-1

Bench: TB

Card Type: GN Testing Instrument: 000017919CE3 (15220)

Bionumber: 0241010103500342

Organism Quantity:

Comments:	

Identification Information	Card: GN	Lot Number: 2410373103	Expires: Dec 4, 2018 12:00 CST			
	Completed: Feb 8, 2018 19:54 CST	Status: Final	Analysis Time: 6.00 hours			
Selected Organism	99% Probability Acinetobacter baumannii					
	Bionumber: 0241010103500342	Confidence: Excellent identification				
SRF Organism						
Analysis Organisms and Tests to Separate:						
Acinetobacter baumannii complex						
Acinetobacter nosocomialis						
Acinetobacter pittii 44C(1),41C(99),						
Acinetobacter baumannii 44C(99),Arg.hydr.(1),41C(99),						
Acinetobacter calcoaceticus Arg.hydr.(99),41C(1),						
Analysis Messages:						
Contraindicating Typical Biopattern(s)						