

Certificate of Analysis / Quality Control

Product: D-Dimer Latex Agglutination Kit

REF: T9515707

Lot:

Exp.:

Results:

QC-ITEM	SPECIFICATION	RESULT
LATEX REAGENT		
Appearance:	white-opaque	
Specificity: 10 Neg. Plasma Controls	no agglutination	
Phosphate Buffer Solution:	no agglutination	
Relative Sensitivity:	0.75 - 1.2	
POSITIVE CONTROL		
Appearance:	clear, pale straw colour	
Titre:	1:8 - 1:32	
HIV-1&2, HBsAg, HCV Test:	negative	
NEGATIVE CONTROL		
Appearance:	clear, pale straw colour	
Titre:	no agglutination	

This lot meets our specifications.

Rev. 2
Date: 2018-07-25