

Certificate of Analysis

Product Name: APTT-XL - Orig:USA
Product Number: 100402
Lot Number: 251074

Expiration Date: 08/31/2014
Storage: 2 - 8°C

Specification	Specification Limit	Result
Appearance	Clear, yellow solution free of particulate matter.	Complies
Vial-to-vial Variability (CV)	$\leq 1.5 \%$	0.7
Level 1 Master Lot*	26 - 34 S	29
Level 1 Master Lot	Lot Number	147251
Level 2 Master Lot*	42 - 56 S	52
Level 2 Master Lot	Lot Number	V39124
Level 3 Master Lot*	64 - 81 S	78
Level 3 Master Lot	Lot Number	V39684
Heparin Level 2*	50 - 64 S	55
Heparin Level 2	Lot Number	V41730
Normal Donors (mean of 3)	Report	25 S
20% Factor VIII / Level 1 Control	≥ 1.10	1.42
20% Factor IX / Level 1 Control	≥ 1.07	1.23
20% Factor XI / Level 1 Control	≥ 1.21	1.41

For in-vitro diagnostic use.

All clot time specifications are for the Sysmex CA1500. They may be different for other photo-optical instruments.

* Control specifications are lot-specific. Before the stated control lots expire or are depleted, new ranges will be established for the replacement lots.

Connie Carroll

Quality Representative

Date: 2012-10-19