

Certificate of Analysis

Product Name: Calcium Chloride, 0.02M - Orig:USA
Product Number: 100304
Lot Number: 785895


Expiration Date: 12/31/2018
Storage: 2 - 8°C

Specification	Specification Limit	Result
Appearance	Clear, colorless solution essentially free of particulate matter	Complies
Performance v. Previous Lot: Difference	$\leq 3 \%$	1

For in-vitro diagnostic use.

Reactive Component: Calcium Chloride

WARNING: Calcium chloride solution contains sodium azide. Sodium azide under acidic conditions yields hydrazoic acid, an extremely toxic compound. Dilute with running water before discarding, and then flush with large volumes of water. These precautions are recommended to avoid deposits in metal piping in which explosive conditions may develop.



Quality Representative

Date: 2016-01-21