



CERTIFICATE OF BATCH ANALYSIS REPORT

Diesse Diagnostica Senese S.p.A hereby Certifies

that the product: **VACU-TEC** Code: **10200/W**

Lot. N.: **8148** Exp. Date: **2019-03**

Prod. Date: **2018-05**

Packed in the pallet N.: **5332**

have been approved after testing the following controls.

ON LINE QUALITY CONTROLS (model CDP not enclosed)

- sodium citrate solution concentration
- sodium citrate volume
- vacuum level (*)
- correct placement of the cap
- no extra body inside the product

QC Passed



(*) Only for VACU-CODE and VACU-TEC tubes

PACKAGING QUALITY CONTROL (model CDP not enclosed)

- control of the number of the product for each pack

Pack N.: **100**

QC Passed



QUALITY CONTROLS ON FINISHED PRODUCT (model CDP not enclosed)

- pack integrity
- data printed on the product label

QC Passed





Production Manager Signature

Date: 2018-05-28



QC Manager Signature

Date: 2018-05-28

Stamp and signature
n.a.

Document Number:

ODV **78571**