

Elisa Certificate of Analysis

This is to certify that the following reagents were tested and approved by Randox Quality Control Department.

Product: Streptomycin Elisa
Catalogue No.: STP3468
Batch No.: 258683
Size: 96 Tests
Kit Expiry: 2014-04
Date of Manufacture: 2012-11

COMPONENT CHECK

Component	Lot No.	Expiry Date	Appearance check
Microtitre Plate	036ST	2014-06	✓
Wash Buffer	4664EA	2014-07	✓
Conjugate	033ST	2014-04	✓
Conjugate Diluent	4678EA	2014-07	✓
One Shot Substrate	128TMB	2016-03	✓
Stop Solution	4703EA	2015-09	✓
Standard 1	041ST	2014-08	✓
Standard 2	042ST	2014-08	✓
Standard 3	043ST	2014-08	✓
Standard 4	044ST	2014-08	✓
Standard 5	045ST	2014-08	✓
Standard 6	046ST	2014-08	✓

ABSORBANCE DATA

STANDARD VALUE(ng/ml)	ABSORBANCE ACHIEVED
0	2.116
0.49	1.639
1.95	1.084
7.81	0.591
31.25	0.274
65.00	0.177

ACCURACY CHECK

CONTROL	TARGET	RANGE	RESULT
QC L	1.95ng/ml	1.37-2.54ng/ml	1.90ng/ml
QC E	7.81ng/ml	5.47-10.15ng/ml	7.61ng/ml

DATE TESTED: 12 November 2012**ANALYST:****QC PASS/FAIL:** Pass**APPROVED BY:**

ABX

