



Certificate of Analysis

Lot Number: 7098

Expiration Date: 2019-01-31

Production Date: 2017-01-26

Product Name: System Wash Reagent

REF: 01-03-0035 (940mL)

Shelf-Life: Two Years (Unopened), 30 Days (Opened)

Manufactured/Distributed By: Trinity Biotech

Storage Requirements: Tightly Closed. Cool, Dry Location. (2-28°C, 36-82°F)

Intended Use: For Trinity Biotech A_{1c} and Hemoglobin Variants IVD Assay Use

Method of Analysis:

Visual Examination

Standard: A clear, colorless solution free from observable particulate matter

Result: A clear, colorless solution is free from observable particulate matter.
Meets specification.

Specific Gravity

Standard: 0.985 – 0.995 Nominal

Result: Meets specification

Approval

Quality Manager:

Date:

1/27/2017