



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Haemophilus influenzae (type c) Catalog Number: 0185 Lot Number: 185-29 Reference Number: ATCC® 9007™* Purity: Pure Passage from Reference: 3	Expiration Date: 2018/12/31 Release Information: Quality Control Technologist: Mary L Bowman Release Date: 2017/2/22
---	---

Performance	
Macroscopic Features: Large, circular, low convex, entire edge, pale tan, glistening Microscopic Features: Gram negative short rods or coccobacilli; may show marked pleomorphism	Medium: Chocolate Method: Gram Stain (1)
ID System: Vitek NH (1) See attached ID System results document.	Other Features/ Challenges: Results V factor (Strips): negative X factor (Strips): negative X V factor (Strips): positive (1) Beta Lactamase (Cefinase Disk): negative H. influenzae antiserum - Type c: positive  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE

Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.



TESTING CERT #2655.01

bioMérieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Feb 17, 2017 09:57 CST
Printed by: mlb

Isolate Group: 185 29-1

Bench: MB

Card Type: NH Testing Instrument: 000017919CE3 (15220)

Bionumber: 5715040647

Organism Quantity:

Comments:	

Identification Information	Card:	NH	Lot Number:	2450169403	Expires:	May 14, 2018 13:00 CDT
	Completed:	Feb 16, 2017 22:31 CST	Status:	Final	Analysis Time:	6.00 hours
Selected Organism	99% Probability Haemophilus influenzae Bionumber: 5715040647 Confidence: Excellent identification					
SRF Organism						
Analysis Organisms and Tests to Separate:						
Analysis Messages: See product information for additional information.						
Contraindicating Typical Biopattern(s)						

Biochemical Details																	
1	ArgA	(+)	2	GGT	-	3	LysA	+	4	dGAL	+	5	LeuA	+	6	ELLM	+
7	PheA	+	8	ProA	-	10	PyrA	-	13	TyrA	+	15	APPA	-	18	dGLU	+
19	GLYG	-	20	dMNE	-	22	dMAL	-	28	SAC	-	33	NAG	-	36	URE	+
39	BGALi	-	40	ODC	-	41	AARA	-	45	PVATE	-	46	PHC	+	47	dMLT	+
51	MTE	-	52	IGLM	-	59	PHOS	+	61	dRIB2	+	62	OPS	+	64	dXYL	+

Installed VITEK 2 Systems Version: 07.01
MIC Interpretation Guideline:
AES Parameter Set Name:

Therapeutic Interpretation Guideline:
AES Parameter Last Modified: