

Certificate of Analysis

Product Name: U.C.R.P. - Orig:USA
Product Numbers: 100301
Lot Number: 335821

Production Date: 10-JAN-14
Expiration Date: 2016-05-31

Specification	Specification Limit	Result
Appearance Cake	Yellow to buff colored solid	Complies
Appearance Reconstituted	Yellow to straw colored, clear to hazy, free of particulates when viewed against dark background	Complies
Solubility	Less than 2 minutes	Complies
Moisture	≤ 3 %	0

For in-vitro diagnostic use.

See directional insert for assay values.

Each unit of source material used in the preparation of this product has been tested by an FDA licensed method and found nonreactive for HBsAg and negative for antibodies to HIV and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS, or other infectious diseases. This product, like all materials of human origin, should be handled as potentially infectious biological material.

Connie Carroll

Quality Representative

Date: 2014-01-20