



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications

Microorganism Name: Streptococcus pneumoniae
Catalog Number: 0947
Lot Number: 947-102**
Reference Number: ATCC® 49619™*
Purity: Pure
Passage from Reference: 3

Expiration Date: 2020/4/30

Release Information:

Quality Control Technologist: Mary L Bowman
Release Date: 2018/5/31

Performance

Macroscopic Features:

Small to medium, circular, translucent, alpha hemolytic colonies.

Microscopic Features:

Gram positive, lancet and occasionally ovoid shaped, cocci typically in pairs, also singly, or short chains.

Medium:

CNA

Method:

Gram Stain (1)

ID System: Vitek GP (1)

See attached ID System results document.

Other Features/ Challenges: Results

- (1) Catalase (3% Hydrogen Peroxide): negative
- Optochin differential: Sensitive (6 mm disks: ≥ 14 mm)
- (1) Oxacillin (1 mcg - Disk Susceptibility): ≤ 12 mm
- (1) Cefotaxime (30 mcg - Disk Susceptibility): 31 - 39 mm
- (1) Ceftriaxone (30 mcg - Disk Susceptibility): 30 - 35 mm

Amanda Kuperus
Quality Control Manager
AUTHORIZED SIGNATURE

****Disclaimer:** The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.



bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed May 18, 2018 07:53 CDT
Printed by: mlb

Isolate Group: 947 102-3

Bench: MB

Card Type: GP Testing Instrument: 000017919CE3 (15220)

Bionumber: 061111324303530

Organism Quantity:

Comments:	

Identification Information	Card: GP	Lot Number: 2420659203	Expires: Sep 16, 2019 13:00 CDT
	Completed: May 17, 2018 17:10 CDT	Status: Final	Analysis Time: 5.00 hours
Selected Organism	97% Probability Streptococcus pneumoniae		
	Bionumber: 061111324303530	Confidence: Excellent identification	
SRF Organism			
Analysis Organisms and Tests to Separate:			
Analysis Messages:			
Contraindicating Typical Biopattern(s)			
Streptococcus pneumoniae SAL(17),			

Biochemical Details

2	AMY	-	4	PIPLC	-	5	dXYL	-	8	ADH1	-	9	BGAL	+	11	AGLU	+
13	APPA	+	14	CDEX	-	15	AspA	-	16	BGAR	+	17	AMAN	-	19	PHOS	-
20	LeuA	+	23	ProA	-	24	BGURr	-	25	AGAL	+	26	PyrA	-	27	BGUR	-
28	AlaA	+	29	TyrA	+	30	dSOR	-	31	URE	-	32	POLYB	+	37	dGAL	-
38	dRIB	-	39	ILATk	-	42	LAC	+	44	NAG	+	45	dMAL	+	46	BACI	-
47	NOVO	-	50	NC6.5	-	52	dMAN	-	53	dMNE	+	54	MBdG	+	56	PUL	-
57	dRAF	+	58	O129R	-	59	SAL	+	60	SAC	+	62	dTRE	+	63	ADH2s	-
64	OPTO	-															

Installed VITEK 2 Systems Version: 07.01

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: