



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Streptococcus equi subsp. zooepidemicus Catalog Number: 0101 Lot Number: 101-51** Reference Number: ATCC® 43079™* Purity: Pure Passage from Reference: 3	Expiration Date: 2019/10/31 Release Information: Quality Control Technologist: Carol J Stanoch Release Date: 2017/12/11
--	--

Performance	
Macroscopic Features: 2 colony types: medium, low convex with flat margin, the other is small, convex; both are circular with entire edge, glistening, translucent, wide zone of beta hemolysis.	Medium: SBAP
Microscopic Features: Gram positive cocci occurring singly, in pairs or short chains	Method: Gram Stain (1)
ID System: Vitek GP (1) See attached ID System results document.	Other Features/ Challenges: Results (1) Catalase (3% Hydrogen Peroxide): negative Streptococcus Antiserum Group C: positive  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE

**Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.



bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Dec 11, 2017 08:47 CST
Printed by: cjs

Isolate Group: 101 51-1

Bench: CS

Card Type: GP Testing Instrument: 000017919CE3 (15220)

Bionumber: 071450365757651

Organism Quantity:

Comments:	

Identification Information	Card:	GP	Lot Number:	2420509103	Expires:	Apr 19, 2019 13:00 CDT
	Completed:	Dec 8, 2017 11:20 CST	Status:	Final	Analysis Time:	3.75 hours
Selected Organism	99% Probability Streptococcus equi ssp zooepidemicus Bionumber: 071450365757651 Confidence: Excellent identification					
SRF Organism						
Analysis Organisms and Tests to Separate:						
Analysis Messages:						
Contraindicating Typical Biopattern(s)						

Biochemical Details																	
2	AMY	-	4	PIPLC	-	5	dXYL	-	8	ADH1	+	9	BGAL	+	11	AGLU	+
13	APPA	+	14	CDEX	-	15	AspA	-	16	BGAR	-	17	AMAN	-	19	PHOS	+
20	LeuA	+	23	ProA	-	24	BGURr	+	25	AGAL	-	26	PyrA	-	27	BGUR	-
28	AlaA	+	29	TyrA	+	30	dSOR	-	31	URE	-	32	POLYB	+	37	dGAL	+
38	dRIB	+	39	ILATk	-	42	LAC	+	44	NAG	+	45	dMAL	+	46	BACI	+
47	NOVO	+	50	NC6.5	-	52	dMAN	+	53	dMNE	+	54	MBdG	+	56	PUL	+
57	dRAF	-	58	O129R	+	59	SAL	+	60	SAC	+	62	dTRE	-	63	ADH2s	+
64	OPTO	+															

Installed VITEK 2 Systems Version: 07.01

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: