



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

|                                                                                                                                                                                                                                    |                                                                                                                                                           |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Specifications</b><br><b>Microorganism Name:</b> Vibrio furnissii<br><b>Catalog Number:</b> 0720<br><b>Lot Number:</b> 720-41<br><b>Reference Number:</b> NCTC 11218<br><b>Purity:</b> Pure<br><b>Passage from Reference:</b> 3 | <b>Expiration Date:</b> 2019/2/28<br><b>Release Information:</b><br><b>Quality Control Technologist:</b> Kieshia L Negen<br><b>Release Date:</b> 2017/5/3 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|

| Performance                                                                                                                                                         |                                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Macroscopic Features:</b><br>Large, circular, convex, entire edge, cream, glistening. Most colonies are opaque but an occasional translucent colony may be seen. | <b>Medium:</b><br>Marine Agar                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Microscopic Features:</b><br>Gram negative straight or curved rod.                                                                                               | <b>Method:</b><br>Gram Stain (1)                                                                                                                                                                                                                                                                                                                                                                  |
| <b>ID System:</b> MALDI-TOF<br>See attached ID System results document.                                                                                             | <b>Other Features/ Challenges: Results</b><br>(1) Oxidase (Kovacs): positive<br>TCBS agar: growth, yellow colonies<br>SBAP (growth from pellet): good growth<br>(1) Purple Broth W/Dextrose: positive with weak production of gas<br><br><br>Amanda Kuperus<br>Quality Control Manager<br>AUTHORIZED SIGNATURE |

Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(1) These tests are accredited to ISO/IEC 17025:2005.



TESTING CERT #2655.01

# Bruker Daltonik MALDI Biotyper Classification Results



## Meaning of Score Values

| Range       | Interpretation                      | Symbols | Color  |
|-------------|-------------------------------------|---------|--------|
| 2.00 - 3.00 | High-confidence identification      | (+++)   | green  |
| 1.70 - 1.99 | Low-confidence identification       | (+)     | yellow |
| 0.00 - 1.69 | No Organism Identification Possible | (-)     | red    |

## Meaning of Consistency Categories (A - C)

| Category | Interpretation                                                                                                                                                                                                                                                                                                                 |
|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (A)      | <b>High consistency:</b> The best match is a high-confidence identification. The second-best match is (1) a high-confidence identification in which the species is identical to the best match, (2) a low-confidence identification in which the species or genus is identical to the best match, or (3) a non-identification. |
| (B)      | <b>Low consistency:</b> The requirements for high consistency are not met. The best match is a high- or low-confidence identification. The second-best match is (1) a high- or low-confidence identification in which the genus is identical to the best match or (2) a non-identification.                                    |
| (C)      | <b>No consistency:</b> The requirements for high or low consistency are not met.                                                                                                                                                                                                                                               |

Sample Name: Vibrio furnissii  
 Sample Description: 0720  
 Sample ID: 720-41  
 Sample Creation Date/Time: 2017-04-07T09:28:45.704KN  
 Applied MSP Library(ies): BDAL, Mycobacteria Library (bead method), Filamentous Fungi Library 1.0, Listeria

| Sample Name    | Sample ID | Organism (best match) | Score Value |
|----------------|-----------|-----------------------|-------------|
| C5<br>(+++)(A) | 720-41    | Vibrio furnissii      | 2.03        |

Comments:

N/A