



Certificate Of Performance

Specimen Description: **Ascaris lumbricoides Parasite Suspension (FP05)**

Lot #: FP05-09-16

Expiration Date: 2020/7/31

Fecal suspension containing *Ascaris lumbricoides* eggs fixed in formalin

Minimum of three *Ascaris lumbricoides* eggs per coverslip area

Eggs are oval, having dimensions of 45-75 μm long by 35-50 μm wide, with a golden to tan, thick, mamillated coat

In addition, *Trichuris trichiura* and *Necator americanus* may also be seen.

Verified By: Tracy A Blenker

Title: Quality Control Technologist

A handwritten signature in black ink that reads "Amanda Kuperus".

Amanda Kuperus
Quality Control Manager
AUTHORIZED SIGNATURE