


Certificate of Analysis

Kit Name: FSH ELISA
Lot # FSS5453

Catalog # FS232F
Expiration Date: 2020-01

Component	Part Number	Lot Number	Expiration Date
Microwells coated with Streptavidin	GS224EP	SA2361	2020-12
FSH Standard 1	FS805SS-1	SS2378	2020-01
FSH Standard 2	FS805SS-2	SS2378	2020-01
FSH Standard 3	FS805SS-3	SS2378	2020-01
FSH Standard 4	FS805SS-4	SS2378	2020-01
FSH Standard 5	FS805SS-5	SS2378	2020-01
FSH Standard 6	FS805SS-6	SS2378	2020-01
FSH Enzyme Conjugate: 1 bottle (ready to use)	FH954ER	ER1167	2020-01
TMB Substrate: 1 bottle (ready to use)	CC956TM	TM0177	2021-05
Stop Solution: 1 bottle (ready to use)	CC604SS	SS0146	2020-12
Wash Concentrate, 20X: 1 bottle	CC601WC	W20111	2020-11

Standard Curve



	OD 450nm	Conc. (mIU/mL)
Standard 1	0.012	0
Standard 2	0.267	5
Standard 3	0.524	10
Standard 4	1.064	25
Standard 5	1.716	50
Standard 6	2.488	100

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

Pass ☐ Reject ☐

Signature: _____

Date: _____

