



Certificate Of Performance

Specimen Description: **Endolimax nana Parasite Suspension (FP12)**

Lot #: FP12-10-2

Expiration Date: 2020/7/31

Fecal suspension containing formalin and Endolimax nana parasite

Minimum of ten Endolimax nana cysts per coverslip

Cysts are oval to spherical with a 5-10 μm diameter, and containing four nuclei. Each nuclei has a karyosome.

Verified By: Carol J Stanoch

Title: Quality Control Technologist

A handwritten signature in black ink that reads "Amanda Kuperus". The signature is written in a cursive, flowing style.

Amanda Kuperus
Quality Control Manager
AUTHORIZED SIGNATURE