

DIESE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI) - Italia - Via delle Rose 10
Tel +39-0577-587111 Fax +39-0577-318690



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc 6
Tel +39-02-4859121 r.a. Fax +39-02-48008530

QUALITY CONTROL CERTIFICATE

KIT NAME: U1-70 RNP - CHORUS

Code 86088	
Lot No. 978	Expiry 2014-07

REAGENTS

Label	Description	Lot.	Exp	No
D3662	U1-70 RNP	978	2014-07	6X6
B3662	CALIBRATORE U1-70 RNP	E13220	2016-05	1
C3662	CONTROL SERUM U1-70 RNP	E13220	2016-05	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	340	≥ 100
Control Serum (AU/mL)	36.2	22.8-60.7

CALIBRATOR

1	0	4	9	7	8	5	4	1	9	3	3	6	0	5	9	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader

YES

☒

Operator Signature

SB



CONTROL SERUM

1	0	4	9	7	8	5	1	2	2	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES

☒

Operator Signature

SB



The test has been carried out in accordance with the package insert

TESTED BY: Sara Benocci

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

SB

DATE 28/10/2013

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.