



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Campylobacter jejuni subsp. jejuni Catalog Number: 0481 Lot Number: 481-96** Reference Number: ATCC® 33291™* Purity: Pure Passage from Reference: 3		Expiration Date: 2019/9/30 Release Information: Quality Control Technologist: Carol J Stanoch Release Date: 2017/12/4
Performance		
Macroscopic Features: Small to medium, circular, convex, entire edge, tan to pinkish, smooth and glistening.; some colonies may be larger and mucoid. Some colonies may exhibit tailing along the streak line.	Medium: Chocolate	
Microscopic Features: Slender, spirally curved, gram negative rods.	Method: Gram Stain (1)	
ID System: MALDI-TOF (1) See attached ID System results document.	Other Features/ Challenges: Results (1) Catalase (3% Hydrogen Peroxide): positive (1) Oxidase (Kovacs): positive Rapid Hippurate Hydrolysis: positive Campy Isolation Agar (42 C, MA): Good growth at 48 hrs. Campy Isolation Agar (35 C, MA): Good growth at 48 hrs. Campy Isolation Agar (25 C, MA): No growth at 96 hrs. Motility (Wet Mount): darting-positive  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE	



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications

Microorganism Name: Campylobacter jejuni subsp. jejuni
Catalog Number: 0481
Lot Number: 481-96**
Reference Number: ATCC® 33291™*
Purity: Pure
Passage from Reference: 3

Expiration Date: 2019/9/30

Release Information:

Quality Control Technologist: Carol J Stanoch

Release Date: 2017/12/4

**Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

(1) These tests are accredited to ISO/IEC 17025:2005.

