



Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

Specifications Microorganism Name: Candida albicans Catalog Number: 0332 Lot Number: 332-109 Reference Number: ATCC® 14053™* Purity: < 0.1% Total Pellet CFU Passage from Reference: 3	Expiration Date: 2019/4/30 Release Information: Quality Control Technologist: Mary L Bowman Release Date: 2017/6/9
---	---

Performance	
Macroscopic Features: Medium to large, white to cream, waxy, circular, convex, entire edge Microscopic Features: Gram positive, spherical, budding yeast cells	Medium: PDA Method: Gram Stain (1)
ID System: Vitek YST (1) See attached ID System results document.	Other Features/ Challenges: Results (1) Germ Tube Test: positive (1) Chlamyospore production: positive  Amanda Kuperus Quality Control Manager AUTHORIZED SIGNATURE

Disclaimer: The last digit(s) of the lot number appearing on the product label and packing slip are merely a packaging event number. The lot number displayed on this certificate is the actual base lot number.

Note for Vitek®: Although the Vitek® panel uses many conventional tests, the unique environment of the card, combined with the short incubation period, may produce results that differ from published results obtained by other methods.

⚠ Refer to the enclosed product insert for instructions, intended use and hazard/safety information.

Individual products are traceable to a recognized culture collection.



(*) The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.



(1) These tests are accredited to ISO/IEC 17025:2005.

TESTING CERT #2655.01

bioMerieux Customer: 1012555
System #: 521711

Laboratory Report

Printed Jun 7, 2017 11:17 CDT
Printed by: mlb

Isolate Group: 332 109-1

Bench: MB

Card Type: YST Testing Instrument: 000017919CE3 (15220)

Bionumber: 6102546065327771

Organism Quantity:

Comments:	

Identification Information		Card: YST	Lot Number: 243397610	Expires: Nov 2, 2017 13:00 CDT
		Completed: Jun 1, 2017 06:18 CDT	Status: Final	Analysis Time: 18.25 hours
Selected Organism		99% Probability Candida albicans		
		Bionumber: 6102546065327771		Confidence: Excellent identification
SRF Organism				
Analysis Organisms and Tests to Separate:				
Analysis Messages:				
Contraindicating Typical Biopattern(s)				

Biochemical Details

3	LysA	-	4	IMLTa	+	5	LeuA	+	7	ARG	+	10	ERYa	-	12	GLYLa	-
13	TyrA	(-)	14	BNAG	-	15	ARBa	-	18	AMYa	-	19	dGALa	+	20	GENa	-
21	dGLUa	+	23	LACa	-	24	MAdGa	+	26	dCELa	-	27	GGT	-	28	dMALa	+
29	dRAFa	-	30	NAGA1	+	32	dMNEa	+	33	dMELa	-	34	dMLZa	-	38	ISBEa	-
39	IRHAa	-	40	XLTa	+	42	dSORa	+	44	SACa	+	45	URE	-	46	AGLU	+
47	dTURa	+	48	dTREa	+	49	NO3a	-	51	IARAa	-	52	dGATa	+	53	ESC	-
54	IGLTa	+	55	dXYLa	+	56	LATa	+	58	ACEa	+	59	CITa	+	60	GRTas	+
61	IPROa	+	62	2KGa	+	63	NAGa	+	64	dGNTa	+						

Installed VITEK 2 Systems Version: 07.01

MIC Interpretation Guideline:

AES Parameter Set Name:

Therapeutic Interpretation Guideline:

AES Parameter Last Modified: