



Certificate Of Performance

Specimen Description: **Taenia species Parasite Suspension (FP04)**

Lot #: FP04-09-15

Expiration Date: 2020/8/31

Fecal suspension containing Taenia eggs fixed in formalin

Minimum of three Taenia eggs per coverslip area

Eggs are round to slightly oval with a 31-43 μm diameter, and contain a 6-hooked embryo or oncosphere inside a thick, striated shell.

In addition, Entamoeba coli cysts may also be seen.

Verified By: Carol J Stanoch

Title: Quality Control Technologist

A handwritten signature in black ink that reads "Amanda Kuperus". The signature is fluid and cursive.

Amanda Kuperus
Quality Control Manager
AUTHORIZED SIGNATURE