

DIESE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI) - Italia - Via delle Rose 10
Tel +39-0577-587111 Fax +39-0577-318690



Pag 1/1
Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

QUALITY CONTROL CERTIFICATE

KIT NAME: CYTOMEGALOVIRUS IgM-CHORUS

Code 81011

Lot No. 004 **Expiry** 2015-04

REAGENTS

Label	Description	Lot.	Exp	No
D2171	Cytomegalovirus IgM	004	2015-04	6X6
B2171	CALIBRATORE Cytomegalovirus IgM	004	2015-04	1
C2171	CONTROL SERUM Cytomegalovirus IgM	004	2015-04	1
G2171	CONFIRMATORY TEST	004	2015-05	5X1


CONTROLS

	Found results	Expected results
Calibrator (mOD)	410	≥ 100
Control Serum (INDEX)	1,5	1,1-1,9
Negative Control (INDEX)	0,0	$<0,7$

CALIBRATOR

0	7	9	0	0	4	7	0	0	0	0	0	7	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature *[Signature]*



CONTROL SERUM

0	7	9	0	0	4	7	1	1	1	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature *[Signature]*



The test has been carried out in accordance with the package insert

TESTED BY: Laura Capitani

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE 17-02-2014

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.