

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI) – Italia – Via delle Rose 10
Tel +39-0577-567111 Fax +39-0577-318890



Pag 1/1
Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6
Tel. +39-02-4859121 r.a. Fax +39-02-48008530

QUALITY CONTROL CERTIFICATE

KIT NAME: CARDIOLIPIN-M - CHORUS

Code	86048		
Lot No.	339	2018-05	2017-05

REAGENTS

Label	Description	Lot.	Exp	No
D3331	CARDIOLIPIN-M	339	2018-05	6X6
B3331	CALIBRATORE CARDIOLIPIN-M	339	2018-05	1
C3331	CONTROL SERUM CARDIOLIPIN-M	339	2018-05	1

CONTROLS

Found results	Expected results
Calibrator (mOD)	440
Control Serum (MPLU/mL)	44,9
	≥ 100
	36,9-100,0

CALIBRATOR

3	3	6	3	3	9	6	4	1	9	3	2	3	2	4	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader

YES



Operator Signature

AS



CONTROL SERUM

3	3	6	3	3	9	6	2	2	4	1	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader

YES



Operator Signature

AS



The test has been carried out in accordance with the package insert

TESTED BY: Sara Benocci

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

Aliaudi

DATE 13/06/2017

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.